

stryker®

Neuro Spine ENT

ENT Surgery



Your Specialty Is Our Priority.

ENT Surgery

Neuro Spine ENT

Saber Drill

5400-120-000

- Tapered attachments, Telescoping burs for line of sight
- Reduce thermal necrosis with integrated irrigation
- 60,000 RPM



CORE Console

5400-050-000




Uni-Directional Footswitch

5100-007-000



TPS Footswitch

5100-008-000



Elite Burs



Zyphr Burs

Attachments

5100-120-490 SD/PD Series Extra Long Plus Straight
 5100-120-482 SD/PD Series Extra Long Angled
 5100-120-480 SD/PD Series Extra Long Straight
 5100-120-472 SD/PD Series Long Angled
 5100-120-470 SD/PD Series Long Straight

5100-120-452 SD/PD Series Medium Angled
 5100-120-450 SD/PD Series Medium Straight
 5100-120-045 Straight Portmann Chuck
 5100-120-044 Angled Portmann Chuck

Accessories

5100-010-272 Long Irrigation Clip
 5100-010-270 Medium Irrigation Clip
 5400-050-001 Irrigation Cassettes 10/Pkg
 5400-121-000 Saber Handswitch

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 1000-927-005 Rev. A
 UnDe/PS.

Copyright © 2007 Stryker
 Printed in USA

Stryker Neuro Spine ENT
 4100 East Milham Avenue
 Kalamazoo, MI 49001 USA
 t: 269 323 7700 f: 269 353 3026
 toll free: 800 253 3210

www.stryker.com

